User's Manual

PRODIGY®

Blood Glucose Monitoring Systems

"Talking"

Model No. 51720

Model No. 51600

***No-Coding
Required!***

Model No. 51810

***"Talking"
No-Coding
Required!***

Model No. 51850

Diagnostic Devices, Inc.

Dear PRODIGY Owner:

Thank you for choosing one of our PRODIGY blood glucose monitoring systems. This manual contains important information about the PRODIGY system you now own. Please read it thoroughly and carefully. A warranty registration card is included with your system. Please complete it and mail it to us.

PRODIGY meters are designed to help you and your health care professionals monitor your blood glucose levels. This manual is designed to help you use the PRODIGY system, and technical support is available twenty four (24) hours a day and seven (7) days a week by calling **1-800-243-2636**. All questions about interpreting the results should be directed to your healthcare professionals.

Caution: Before using any product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

This manual is written for four (4) models of the PRODIGY meter: PRODIGY Advance, PRODIGY Audio, and PRODIGY Autocode talking meter and non-talking meter.

All four (4) PRODIGY models have the latest technology for blood glucose monitoring, are easy to use, and will give fast and accurate results with a minimal sample of blood. All four (4) PRODIGY meters have a large screen display for easy reading, and are small and lightweight for portable convenience. All the PRODIGY meters allow Alternative Site Testing (AST), as well as have memory and data management capabilities, with free software available at www.prodigymeter.com, to give you and your healthcare professionals powerful graphic tools to manage your diabetes.

PRODIGY Audio and PRODIGY Autocode Talking also include a speaking (talking) function, which is an acoustic aid for using the meter and hearing test results. The PRODIGY Autocode meters (talking and non-talking) are no-coding meters, allowing you to save time and avoid human error due to improper coding.

IMPORTANT SAFETY INSTRUCTIONS READ THIS BEFORE USING

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
2. Use the device only for the intended use described in this manual.
3. Do not use strips and control solutions which are not supplied by the manufacturer.
4. Do not use the device if it is not working properly, or if it has suffered any damage.
5. Before using any product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

KEEP THESE INSTRUCTIONS

Table of Contents

Important Safety Instructions	3
Important Information	5
About Alternative Site Testing (AST)	5
Introduction of the System	7
Intended Use	7
Principle of Measurement	7
Contents of the System	8
Appearance and Key Function of the Meter	10
Meter Display Segments	11
Speaking Function	12
Appearance of the Test Strip	13
Important Test Strip Information	14
Before Testing	15
Checking the Display	15
Calibration	15
About PRODIGY Control Solution	17
Important Control Solution Test Information	17
How to Do a Control Solution Test	18
Testing Your Blood	21
Getting a Drop of Blood	21
Step-by-Step Test Procedure	24
Using the Meter Memory	26
Viewing Results on a Personal Computer	28
Comparing Meter and Laboratory Results	29
Setting the Meter and Deleting the Memory	30
Caring for Your Meter	35
Maintenance	35
Battery Replacement	35
Displaying Messages and Problem-Solving Guide	37
Special Messages	38
Error Messages	39
Troubleshooting Guide	42
Specifications	44
Expected Test Results	44
Symbols Information	47

Important Information

Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.

Elevated blood triglycerides, reducing substances such as uric acid and ascorbic acid at normal blood concentration, or acetaminophen, dopa, methyldopa, L-dopa, and tolbutamide at normal blood concentrations should not significantly affect the test results.

If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.

Test results below 60 mg/dL (3.3 mmol/L) indicate low blood glucose (hypoglycemia). Test results greater than 240 mg/dL (13.3 mmol/L) indicate high blood glucose (hyperglycemia). If you get results below 60 mg/dL or above 240 mg/dL, repeat the test, and if the results are still below 60 mg/dL (3.3 mmol/L) or above 240 mg/dL (13.3 mmol/L), consult your healthcare professional immediately.

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.

Please refer to your test strip package insert for additional important information.

About Alternative Site Testing (AST)

Important: There are important limitations to AST. Please consult your healthcare professionals before you do AST.

What is AST?

Alternative Site Testing (AST) means you can use parts of the body other than your fingertips to check your blood glucose levels. The PRODIGY blood glucose meters allow you to test on your palm, forearm, upper arm, calf, or thigh.

What is the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, nerve endings are not so numerous, and you will not feel as much pain as you will experience at the fingertip.

When to use AST?

Food, medication, illness, stress, and exercise all can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Therefore, if you are testing your blood glucose level during or immediately after a meal, physical exercise, or stressful event, *take the blood sample from your fingertip only.*

Use AST only:

- Two hours or more after your last meal.
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You have reason to believe you have hypoglycemia or hyperglycemia.
- Your routine glucose results are often fluctuating.
- You are pregnant.

To increase the accuracy when using AST, rub the puncture site before extracting blood.


Figure 1

Introduction of the System

• **Intended Use**

The system is intended for use outside the body (in vitro diagnostic use only). It should be used only for testing glucose (sugar) and only with fresh capillary whole blood samples (from the fingertip, the palm, the forearm, the upper arm, the calf and the thigh). The system is intended for use in the home and in clinical settings. It should not be used for the diagnosis of diabetes or for the testing of newborns.

Alternative Site Testing (AST) in this system can be used only during steady-state blood glucose conditions described in the “About AST” section.

• **Principle of Measurement**

The test is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures electric current and displays the corresponding blood glucose level. The strength of the electric current produced by the reaction depends on the amount of glucose in the blood sample.

Contents of the System

This Manual is written for ALL four (4) PRODIGY systems: PRODIGY Advance, PRODIGY Audio, PRODIGY Autocode Talking and Non-Talking meters. Each model is available either as a meter alone or as a kit. Each kit includes the supplies needed to start testing your blood glucose, i.e., test strips, control solution, lancets and lancing device. If you have a meter only, you can purchase these supplies from your provider. Please check the “REF” number marked on the side panel of the box to see if you have purchased a “Meter Kit” or a “Meter”. Please review the contents of your purchase to confirm that all of the components are included as listed below:

Your “Meter Kit” includes:

- Meter with 2-AAA Batteries
- Carrying Case
- 1- 4ml Bottle of Control Solution
- 10 - Prodigy Test Strips
- 10 - Sterile Lancets
- Lancing Device (with clear cap)
- Complete Instructions:
 - Manual
 - Log Book
 - Quick Reference
 - Warranty Card

Your “Meter” includes:

- Meter with 2-AAA Batteries
- Carrying Case
- Complete Instructions:
 - Manual
 - Log Book
 - Quick Reference
 - Warranty Card

The PRODIGY meters, test strips, and control solutions have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Use only PRODIGY control solutions, and use only the proper PRODIGY test strips for your PRODIGY meter model. PRODIGY Advance and PRODIGY Audio models use PRODIGY test strips REF. 52400 and REF. 52100. Both PRODIGY Autocode meters use PRODIGY Autocode test strips REF. 52800 and REF. 52810.

Note: Check your system to be sure that it is unopened prior to use and that it contains all parts listed. If your product has been opened before use or does not contain all parts listed above, please return your system to the place of purchase.

Appearance and Key Function of the Meter

Test Slot

Insert the test strip here. The meter will turn on automatically.

Talking Symbol

Appears on PRODIGY Audio and PRODIGY Autocode to confirm audio function.

LCD Display

Guides you through the test using symbols and simple messages.

Main Button

Press to turn the meter on or perform other functions described in this manual.

Code Button

On PRODIGY Advance and PRODIGY Audio, used to code meter. PRODIGY Autocode meters do this automatically.

Data Port

For Interface Cable connection to your computer.

Set Button

(in the back) Located in the battery compartment, used to set up the meter.


Figure 2


Figure 3


Figure 4

Meter Display Segments

Test Strip Symbol

Appears when the meter is turned on.

Blood Drop Symbol

Flashes when it is ready to apply the sample.

Code

Shows current code number for PRODIGY Advance and PRODIGY Audio. PRODIGY Autocode meters do this automatically.

Temperature Symbol

Appears when ambient temperature is outside of operating range.

Date

Time

CTL Symbol

Appears when doing a control test. Result will not be stored in memory.

Day Average

Indicates that the displayed result is an average.

Memory Symbol

Appears when you review the memory.

Voice Symbol

Shows if audio function is on, for PRODIGY Audio and PRODIGY Autocode (Talking).

Test Result Area

Displays glucose results. A decimal point appears when the measurement unit is mmol/L.

Low Battery Symbol

Appears when the battery power is low.

Face and Low/High Symbol

Indicates that test result is outside of reference range.

Measurement Unit

Appears with the test result either in mg/dL or mmol/L.

Ketone Warning

Appears when the test result is equal or higher than 240mg/dL (13.3mmol/L).


Figure 5

Speaking “Talking” Function (Model No. 51720, 51850 only)

PRODIGY Audio and PRODIGY Autocode Talking meters will “speak” or “talk” step-by-step through your test. The following table shows you when the meter “speaks” or “talks” and what it says.

WHEN does the meter speak?	WHAT does the meter say?
When the meter is turned on.	<i>Welcome music plays.</i> <i>Thank you for using PRODIGY.</i>
When room temperature is detected.	<i>The room temperature is (number) degree Fahrenheit/ degree Celsius.</i>
When room temperature is outside operating range, which is 50°F~104°F (10°C-40°C).	<i>Room temperature out of range, unable to measure.</i>
After inserting test strip and completing code selection.	<i>The code number is (number).</i> *Only for PRODIGY Audio.
When the meter is ready to test. (💧 symbol appears on display)	<i>Please apply blood into the strip.</i>
When the test is completed. (The result appears on display)	<i>The blood glucose is (number) milligrams per deciliter/ millimoles per Liter</i>
When the test result is outside measurement range, which is 20-600 mg/dL.	<i>The blood glucose is out of range, unable to measure.</i>
When you turn off the meter.	<i>Music plays.</i>

Appearance of the Test Strip

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent channel of the test strip and is automatically drawn into the reaction cell.

The test strip consists of the following parts:

Absorbent Channel

Apply a drop of blood here. The blood will be drawn automatically.

Confirmation Window

This is where you confirm whether enough blood has been applied to the strip absorbent channel.

Test Strip Handle

Hold this part to insert the test strip into the meter slot.

Contact Bars

Insert this end of the strip into the meter. Push it in firmly until it will go no further.

Figure 6

Please see pages 21~25, “Testing Your Blood”, for complete instructions.

Important Test Strip Information

Each model of the PRODIGY meters uses its own PRODIGY test strips. The PRODIGY Advance and the PRODIGY Audio meter will only work with PRODIGY test strips REF. 52400 and REF. 52100. The PRODIGY Autocode meters will only work with PRODIGY test strips REF. 52800 and REF. 52810.

- Store test strip packages in a cool, dry place between 39.2°F~104°F (4°C~40°C). Keep away from direct sunlight and heat. Do not refrigerate.
- Store your test strips in their original vial only. Do not transfer them to a new vial or any other container.
- After removing a test strip from the vial, immediately replace the vial cap and close it tightly.
- With clean, dry hands, you may touch the test strip anywhere on its surface when removing it from the vial or inserting it into the meter.
- Use each test strip immediately after removing it from the vial.
- Apply only a blood sample or control solution to the strip absorbent channel. Applying other substances to the strip absorbent channel will cause inaccurate results.
- Write the discard date on the vial label when you first open it. Discard remaining test strips 90 days after first opening date.
- Do not use test strips beyond the expiration date printed on the package.
- Do not bend, cut, or alter a test strip in anyway.

Warning: Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard. Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.

Before Testing

Checking the Display

Each time you insert a test strip the meter displays “CH” and “

For PRODIGY Audio only, you will hear this voice message: **“The code number is (number).”**

Calibration

Note: PRODIGY Autocode meters do not require calibration. If you are using either of the PRODIGY Autocode meters, skip this section.

If you are using PRODIGY Advance or PRODIGY Audio, you will need to calibrate the meter when you first use it or when you begin using a new vial of test strips.

To calibrate the Meter, press the Code Button, and follow these Steps:

A) To Calibrate When the Meter is Off, Press the Code Button.

The screen displays “Code” and “SEL”, and then a flashing code number. Press the **Code** Button and repeat until the code number on the Meter matches the Code Number printed on the test strip vial. Then press the **M** button to set the code. “☺” and “OFF” then appear, which indicate that calibration is completed and the meter is turned off.


Figure 7

Figure 8



Figure 9

B) To Calibrate When the Meter Has a Test Strip Inserted

Step 1 – Insert the Test Strip

Start with the meter off. Insert a test strip into the test slot. LCD displays “CH” and “” first, the ambient temperature next, and finally a flashing code number.

Step 2 – Match the Code Number

Compare the flashing number with the code number on the test strip vial. If the two numbers match, press the **M** button or wait 5 seconds to fix the code. Then “” appears, indicating that you can directly begin testing. If they do not match, proceed to Step 3.

Step 3 – Select a Correct Code

While the code number is still flashing, press the **Code** button to advance one number. When the correct code appears, press the **M** button or wait 5 seconds to fix the code. Then “” appears. You can begin testing. If you would like to turn off the meter, remove the test strip or wait 2 minutes for the meter to automatically shut down.

Caution: If the code number displayed on the meter does not match the number printed on the vial, test results may be inaccurate.

Figure 10


Figure 11


Figure 12

About PRODIGY Control Solution

PRODIGY® Control Solutions contain a known amount of glucose, which reacts with the test strips. By testing your Control Solution and comparing the results with the expected range printed on the test strip vial label, you can make sure that the meter and the test strips are working properly together. It is very important that you do this simple check routinely to make sure you get accurate results.

How often should the control solution test be performed?

- When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- To continue to get accurate results, you should perform a test with control solution at least once a week.
- Also perform a control solution test whenever you begin using a new vial of test strips, when your test strips are exposed to extreme environmental conditions (See “*Important Test Strip Information*” section of this manual), if you drop the meter, or if you change the batteries.

Important Control Solution Test Information

- Use only PRODIGY control solutions.
- Check the expiration date on the control solution bottle. Do not use if expired.
- Control solution, meter, and test strips should come to room temperature (68-77°F / 20-25°C) before testing.
- Shake the bottle well before use. Discard the first drop of control solution. After squeezing, wipe off the dispenser tip to avoid contamination. These steps ensure you to get a good sample and an accurate result.

- Use within a period of 90 days from the date that you first open it. Record the discard date (date opened plus 90 days) on the control solution bottle. Discard after 90 days.
- Store the control solution tightly closed at temperatures below 86°F (30°C). Do not refrigerate.

Note: The control solution range printed on the test strip vial is for PRODIGY control solution only. It is used to test meter and test strip performance. It is not a recommended range for your blood glucose level.

How to Do a Control Solution Test

Step 1 – Insert Test Strip

Insert a test strip with contact bars end first into the test slot. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically.

After the “” symbol appears on display, press the **M** button, and “**CTL**” appears on the display. With the “**CTL**” symbol on the display, the meter will not store your test result in memory. If you decide not to perform a control solution test, press the **M** button again and the “**CTL**” symbol will disappear.

Note: Every time you perform a control solution test, you should enter into the “**CTL**” test mode so that the test result will not be stored in the meter memory.

Figure 13

Figure 14

Step 2 – Apply Control Solution

Shake the control solution bottle well.

Remove the cap. Squeeze the bottle, discard the first drop, and wipe off the dispenser tip with a clean tissue paper or cotton. Squeeze the bottle again to get a second drop onto a clean non-absorbent surface or on your fingertip first. Then apply the drop to the test strip absorbent channel, (where it meets the narrow channel) until the drop is drawn into the test strip and the confirmation window is filled. The meter begins to count down.




Figure 15

Caution: To avoid contaminating the control solution with the content of the test strip, you have to place a drop of control solution on a clean surface or on your fingertip first. Then touch the strip to the drop. **DO NOT DIRECTLY APPLY CONTROL SOLUTION ONTO A STRIP.**

Correct ✓ – Figure 16

Incorrect – Figure 17

Step 3 – A Result Appears in 6 Seconds

After counting to 0, the control solution test result appears. Compare the result with the range printed on the test strip vial. The result should fall within the printed range.

For PRODIGY Audio and PRODIGY Autocode Talking only. You will hear a voice message, ***“The blood glucose is (number) (unit)”***. (Ex: 137 mg/dL)

Out-of-Range Results

If test results fall outside the range printed on the test strip vial, check the “Troubleshooting Guide” located in the “Displaying Messages and Problem-Solving Guide” section and repeat the test.

If you continue to get out-of-range results, it means that the system may not be working properly. DO NOT use the system to test your blood glucose level. If you are unable to resolve the problem, contact Technical Support at **1-800-243-2636**.


Figure 18

Testing Your Blood

Getting a Drop of Blood

Step 1 – Insert a Lancet in the Adjustable Lancing Device

Twist the Lancing Device cap to remove it. Insert a lancet into the lancet holder and push down firmly until it is fully seated. Do not twist the lancet.

Remove the protective cap from the lancet.

Warning: To reduce the chance of infection:

- Never share a lancet or the Lancing Device with anyone.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the Lancing Device.

Replace and screw the lancing device cap until it is snug but not too tight.

The adjustable tip offers 5 levels of skin penetration. To select the desired depth, twist the adjustable tip in either direction until the number lines up with the arrow. To select the best depth: 1-2 for soft or thin skin, 3 for average skin, 4-5 for thick or calloused skin.


Figure 19


Figure 20


Figure 21

•AST–Alternative Site Testing

A clear cap is available when you purchase a PRODIGY Meter Kit or a PRODIGY Lancing Device, it makes it easier for you to get a drop of blood for AST. When you want to obtain blood from sites other than the fingertip, replace the clear cap instead of the lancing device cap. Screw the clear cap until it is snug but not too tight, and then go to step 2.

Step 2 – Cock the Lancing Device

Slide the ejection/cocking control back until it clicks. If it does not click, the Lancing Device may have been cocked when the lancet was inserted.

Step 3 – Wash Your Hands and the Puncture Site

Use warm, soapy water. Rinse and dry thoroughly.

Step 4 – Select and Lance a Puncture Site

•Fingertip

Hold the Lancing Device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.


Figure 22


Figure 23

• Sites Other Than Your Fingertip

Please refer to the “About AST” section. Please consult your health professional before obtaining blood from sites other than your fingertip.

Warning:

- The system requires a very small blood drop to perform a test. You can obtain it from a fingertip. Choose a different spot each time you test. Repeated punctures in the same spot may cause soreness and calluses.
- Before you decide to do palm testing, please consult your health professional.


Step 5 – Gently Massage the Area

Do not smear the blood sample. To obtain best accurate results, the first drop of blood should be discarded with a clean tissue paper or cotton. Use the second drop of blood for testing. The required blood sample is 0.6 microliter in volume (● actual size).


Step 6 – Remove the Lancet

Take the lancet out carefully. Place the protective cap back on the lancet exposed tip. Always use caution when removing the lancet. Discard the lancet according to your local regulation.

Warning: The first drop of blood usually contains tissue fluid and serum, which may affect the test result. It should be discarded.

Correct – Figure 24

Incorrect – Figure 25

Step-by-Step Test Procedure

Step 1 – Insert the Test Strip

Insert a test strip with contact bars end first into the test slot. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically.

Step 2 – Apply Sample

For PRODIGY Audio and PRODIGY Autocode Talking only. You will hear a voice message, **“Please apply blood into the strip”**.

When the meter shows “”, obtain a drop of blood using the adjustable Lancing Device. The blood sample must be at least 0.6 microliter in volume (● actual size) or you may get an inaccurate test result. Apply blood to the strip absorbent channel where it meets the narrow channel. Do not push your finger against the test strip or try to apply a smeared sample. The Test Strip confirmation window should be completely filled before the meter begins to count down. If not, do not try to add more blood to the test strip. Discard the test strip and retest with a new one. If you have trouble filling the test strip, please contact Technical Support for assistance.




Figure 26

Correct – Figure 27a

Incorrect – Figure 27b


Figure 28

Step 3 – Obtain an Accurate Result in 6 Seconds

Your blood glucose test result appears after the meter counts down to 0. On the PRODIGY Audio and PRODIGY Autocode talking meters, the speaker announces the test result. This blood glucose result is automatically stored in the meter memory. Turn the meter off by removing the test strip. Discard the used test strip carefully to avoid contamination.

For PRODIGY Audio and PRODIGY Autocode Talking only. You will hear a voice message, **“The blood glucose is (number) (unit)”**. (Ex. 70 mg/dL).

Note: If you do not apply a blood sample within 3 minutes, the meter will automatically turn off. You must remove the test strip and re-insert it again to restart the test procedure.

Figure 29

Using the Meter Memory

Your Meter stores the 450 most recent blood glucose test results with date and time in its memory. It also provides you with 7, 14, 21, 28, 60, and 90 day averages of your blood glucose test results. You can review the test results in memory with these easy steps.

Step 1 – Enter the Memory Mode

With the meter turned off, press the **M** button twice. The 7-day average will appear, indicating that you are in the memory mode. If you continue to press the **M** button, the 14, 21, 28, 60, and 90 day averages will appear in order. You can then review the last 450 individual test results in memory.

When using the meter for the first time, “---” appears, showing that there are no test results in memory.

The 7 day average is calculated from the blood glucose results obtained during the last 7 days.

It also indicates how many blood glucose tests have been performed within this period, e.g., 15 (15 tests in the last 7 days).

The 14 day average is calculated from the blood glucose results obtained during the last 14 days. It, too, indicates how many blood glucose tests have been performed, e.g., 27 (27 tests in the last 14 days). The 21, 28, 60 and 90 day averages show the same information.


Figure 30


Figure 31

Step 2 – Recall Test Results

After the 90 day average, the most recent test result with date and time will be shown. Press the **M** button once and the next most recent test result will appear. Each time you press and release the **M** button, the meter will recall up to your last 450 test results in order. When the memory is full, the oldest result is dropped as the newest is added.

When reaching the last set of results, **“End”** will be shown on the LCD display, and the meter will be turned off automatically.

Step 3 – Exit The Memory Mode

Press and hold the **M** button for 3 seconds to turn off the meter.

Note: If no button is pressed for 2 minutes, the meter will show **“OFF”** and turn off automatically.


Figure 32

Figure 33

Viewing Results on a Personal Computer

Results in memory can be transmitted to your personal computer. Prodigy Diabetes Management System Software and an Interface Cable are required before installation. The software can be downloaded free from www.prodigymeter.com. The Interface Cable is an optional accessory. To learn more about Prodigy Diabetes Management System Software or to obtain an Interface Cable (sold separately), please contact your local diabetes supplies service.

Step 1 – Install Software

Install Prodigy Diabetes Management System Software on your Computer by following the instructions provided on our website:

Go to <http://www.prodigymeter.com> and click on "Data Management".

Step 2 – Connect to a Personal Computer

Connect the Interface Cable to your computer. With the meter turned off, connect the Interface Cable to the Data Port of the meter. "PC" will appear on the display, indicating that the meter is ready to transmit data.

Step 3 – Transmit Data

Follow the instructions provided in the software to transmit data. Results transmitted will include date and time. Remove the Interface Cable and the meter will automatically turn off.


Figure 34

Note: While the meter is connected to the PC, it is unable to perform a blood glucose test.

Comparing Meter and Laboratory Results

The result you obtain from your meter may differ somewhat from your laboratory results due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. (See test strip package insert for typical accuracy and precision data, and for important information on limitations.) To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before you go to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- It is strongly recommended to fast for at least eight (8) hours before doing comparison tests.
- Take your meter with you to the lab.

While at the lab:

- Make sure that the samples for both tests (the meter test and the lab test) are taken and tested within 15 minutes of each other.
- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a graytop test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication, or experienced stress¹. In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test². Therefore, it is best to fast for eight (8) hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a meter result to be different from a laboratory result.

References: 1) Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

2) Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia:W.B. Saunders Company (1994),959.

Setting the Meter and Deleting the Memory

Your meter comes with the time, date, unit of measurement and unit of temperature preset. However if you need to change the time, or if you replace the batteries, you may need to enter the set-up mode.

To set the time, you must first enter the set-up mode.

Start with the meter off. Then press the **Set** button located in the battery compartment. The meter is now in the set-up mode.

Step 1 – Set the Year

The year appears first, with the number flashing. Press and release the **M** button to advance one year. To move faster, hold the **M** button down. With the correct year on the display, press the **Set** button and the date will appear on the display with the month segment flashing.

Step 2 – Set the Month

Press and release the **M** button until the correct month appears. To move faster, hold the **M** button down. With the correct month on the display, press the **Set** button and the day segment will start flashing.

Figure 35


Figure 36

Step 3 – Set the Day

Press and release the **M** button until the correct day appears. To move faster, hold the **M** button down. With the correct day on the display, press the **Set** button and the time will appear on the display with the hour segment flashing.

Step 4 – Set the Hour

Press and release the **M** button to advance one (1) hour. To move faster, hold the **M** button down. With the correct hour on the display, press the **Set** button and the minutes setting will start flashing.


Figure 37


Figure 38

Step 5 – Set the Minutes

Press and release the **M** button to advance one minute. To move faster, hold the **M** button down. With the correct minute on the display, press the **Set** button and the current unit of measurement will start flashing.

Step 6 – Select mg/dL or mmol/L

Press and release the **M** button until the unit of measurement you are choosing appears on the display. Press the **Set** button and the current unit of temperature will start flashing.

Note: Day averages are calculated from results obtained during the 7, 14, 21, 28, 60 and 90 days preceding the current date and time settings. When the date and time are changed, the 7, 14, 21, 28, 60 and 90 day averages may change.

Your meter can display test results in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L). The mg/dL unit is standard in the United States. The mmol/L unit is standard in Canada.

Step 7 – Select °F or °C

Press and release the **M** button until the unit of temperature you are choosing appears on the display. Press the **Set** button and the meter will display “dEL” with flashing “M” symbol.

Figure 39


Figure 40


Figure 41

Step 8 – Delete Memory

When the “dEL” and flashing “M” symbol appear on the display, you can choose to clear the memory in the meter. If you do not want to clear the memory, press the **Set** button again to skip this step. If you want to clear ALL memory, press the **M** button, and both “dEL” and “M” will flash. Press the **M** button again to delete ALL memory. The meter then displays “☺”, which means that the memory is deleted.

If you have a PRODIGY Audio or PRODIGY AutoCode Talking Meter, please press the **Set** Button to proceed to the next step. If not, you have completed set-up. Press the **M** or **Set** button to turn off the meter. “OFF” is displayed before shut down.

Step 9 – Select Speaking “Talking” Function

If you have a PRODIGY Audio or PRODIGY AutoCode Talking Meter, the meter displays “VOL”, “🔊” and a flashing number. Press the **M** button to select the speaking volume from 0 to 7. Then press the **Set** button to proceed to the next step.

Number 0 indicates that the speaking function is turned off, where “🔊” does not display on LCD during testing. Numbers 1 to 7 indicate speaking volume from low to high, where “🔊” displays on LCD during testing.


Figure 42



Figure 43

Figure 44

Step 10 – Select Language

If you have a PRODIGY Audio or PRODIGY Autocode Talking meter, press the **Set** button once again to select your language. “**L1**” or “**L2**” and “” will appear. Press the **M** button to select English which is “**L1**” by default or the alternate language which is “**L2**”.

Set-up is now completed. Press the **Set** button again to turn off the meter. “**OFF**” is displayed before shut down.

Figure 45

Figure 46

Caring for Your Meter

Maintenance

Your meter does not require special maintenance. If no blood or control solution comes in contact with the meter, there is no special cleaning required. Take care to avoid getting dirt, dust, blood, control solution or water inside the meter through the test port or data port. Store the meter in its carrying case after each use. A cloth dampened with water and mild detergent can be used to wipe off the outside of the meter. Your meter is a precision instrument. Please handle it with care.

Battery Replacement

Your meter comes with two (2) 1.5V AAA alkaline batteries. The meter will alert you when the power is getting low by displaying two different messages:

1. When the “” symbol appears alone on the display, the meter is functional and the result remains accurate, but you should change the battery as soon as possible.
2. When the “” symbol appears together with low and E-b symbols on the displays, the batteries do not have enough power for a test. You must change the batteries before using the meter.

Figure 47


Figure 48

To replace the batteries, make sure that the meter is turned off.

STEP 1 – Press the buckle on the battery cover and lift up to remove the cover.

Step 2 – Remove the old batteries and replace with two (2) 1.5V AAA new alkaline batteries.

Step 3 – Close the battery cover.

Note: Replacing the batteries does not affect the meter’s memory (previous test results stored in memory). However, the time, date and units settings may need to be updated.

As with all small objects, the batteries should be kept away from small children. If the batteries are swallowed, seek medical assistance immediately.


Figure 49


Figure 50


Figure 51

Display Messages and Problem-Solving Guide

The following is a summary of some display messages. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under the “Action” section.

If you follow the actions recommended but the problem is not resolved, please call Technical Support at 1-800-243-2636 for assistance.

Special Messages


Special symbols and messages that appear together with your test result.


Message	What it means
 <p>Figure 52</p>	<p>“☺” appears when your result is within the reference range from 70 to 120 mg/dL (3.9 and 6.6 mmol/L).</p>
 <p>Figure 53</p>	<p>“low ☺” appears when your result is between 20 and 69 mg/dL (1.1 and 3.8 mmol/L). This indicates that the result is below the reference range.</p>
 <p>Figure 54</p>	<p>“Lo” appears when your result is below the measurement limit, which is less than 20 mg/dL (1.1 mmol/L).</p>

“Lo” or “low ☺” symbol indicates hypoglycemia (low blood glucose.) You should immediately consult your healthcare professional.

Message	What it means
 <p>Figure 55</p>	<p>“high ☹️” appears when your result is equal to or greater than 120 mg/dL (6.6 mmol/L). This indicates that the result is higher than the reference range.</p>
 <p>Figure 56</p>	<p>“KETONE?” appears together with “high ☹️”, if your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates there is a possibility of ketone accumulation if you are a Type 1 diabetic. You should immediately consult your healthcare professional.</p>
 <p>Figure 57</p>	<p>“Hi” appears when your result is above the measurement limit, which is higher than 600 mg/dL (33.3 mmol/L). You should immediately consult your healthcare professional.</p>

Error Messages

Message	What it means	Action
 <p>Figure 55</p>	<p>Appears when the batteries cannot provide enough power for a test.</p>	<p>Replace the batteries immediately.</p>
 <p>Figure 56</p>	<p>Appears when inserting a used test strip.</p>	<p>Test with a new test strip. If the problem persists, please call Technical Support at 1-800-243-2636.</p>
 <p>Figure 57</p>	<p>Appears when the environmental temperature is below the system operation range.</p>	<p>System operation range is 50°F~104°F (10°C~40°C). Repeat the test after the meter and test strip is within the above temperature range.</p>
 <p>Figure 58</p>	<p>Appears when environmental temperature is above the system operation range.</p>	

Message	What it means	Action
 <p>Figure 59</p>	<p>The strip has been removed after applying blood to the absorbent channel.</p>	<p>Re-test with a new test strip.</p>
 <p>Figure 60</p>	<p>Problem with calibration.</p>	<p>Check if the code number on the meter display matches the code number on the test strip vial and re-test. If the problem persists, please call Technical Support at 1-800-243-2636.</p>
 <p>Figure 61</p>	<p>Problem with the meter.</p>	<p>Review the instructions and re-test with a new test strip. If the problem persists, please call Technical Support at 1-800-243-2636.</p>

Troubleshooting Guide

The meter does not display a message after inserting a test strip.

Probable Cause	What to Do
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end first and facing up.
Defective meter.	Please call Technical Support for assistance at 1-800-243-2636.

The test does not start after applying the sample.

Probable Cause	What to Do
Insufficient blood sample.	Repeat the test using a new test strip with a larger blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shut-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when the “  ” symbol appears on the display.
Defective meter.	Please call Technical Support for assistance at 1-800-243-2636.

Troubleshooting Guide Cont.

If the Control Solution test result is out of range.

Probable Cause	What to Do
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Improper code number.	Check if the code number on the meter display matches the code number on the test strip vial.
Control Solution bottle not shaken well.	Shake the Control Solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date or the discard date of the control solution.
Control Solution that is too warm or too cold.	Control Solution, meter, and test strips should come to room temperature (68-77°F/20-25°C) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Meter malfunction.	Please call Technical Support for assistance at 1-800-243-2636.

Specifications

Dimensions & Weight: 3.79in.(L) x 1.79in.(W) x 1in.(H) 2.76oz.
96mm(L) x 45mm(W) x 23mm(H), 78.19g

Power source: Two (2) 1.5V AAA alkaline batteries

Battery life: Over 1,000 tests

Display: LCD

Memory: 450 measurement results with date and time

External output: RS232 PC interface

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Auto turn-off after 3 minutes without action

Temperature warning

Operating condition: 50°F~104°F (10°C~40°C), below 85% R.H. (non-condensing)

Storage/Transportation condition: 39.2°F~104°F (4°C~40°C), below 85% R.H.

Measurement Units: Either mg/dL or mmol/L

Measurement Range: 20~600mg/dL (1.1~33.3mmol/L)

The specifications may be changed without prior notice.

Performance Characteristics

- **Accuracy:** ± 15 mg/dL when glucose < 75 mg/dL
 $\pm 20\%$ when glucose ≥ 75 mg/dL
 - **Precision:** This study shows the CV (correlation variation) is less than 5%.
 - **The device has certified to meet the following standards:**
98/79/EC, IEC 60601-1, IEC 61010-1, IEC 60601-1-2,
IEC61326, and ISO 15197
-


Expected Test Result

Time of day	Plasma glucose range for people without diabetes (mg/dL)/(mmol/L)	Your target range (mg/dL)/(mmol/L)
Fasting and before meal	Less than 110 / 6.1	_____ (mg/dL)/(mmol/L)
2 hours after meals	Less than 140 / 7.8	_____ (mg/dL)/(mmol/L)
Bedtime	Not specified	_____ (mg/dL)/(mmol/L)
Between 2 AM and 4AM	Not specified	_____ (mg/dL)/(mmol/L)

Source: ADA Clinical Practice Recommendations 2003
Please work with your doctor to determine a target range that works best for you.

The PRODIGY® Blood Glucose Monitoring Systems are designed to help you and your healthcare professionals manage your diabetes. You must always rely on your healthcare professionals to interpret your test results and to decide how to treat your diabetes.

Symbols Information

Symbol	Referent
	Do not re-use. Single use only.
	Consult operating instruction
	Keep away from sunlight
	Keep dry
	Temperature limitation
	Use by
	Date of manufacture
LOT	Batch code
REF	Catalogue number
SN	Serial number
CONTROL	Control

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Charlotte, NC 28269 U.S.A.