

## Using your ACCU-CHEK® Multiclix Lancet Device

### Testing From Your Fingertip

**1. Remove the blue cap.**

Note the position of the thumb and index finger in the picture.



**2. Insert a new lancet drum**

blue ring first, until it stops and clicks into place.

*Note: Before inserting a lancet drum, be sure the device is not primed. (The release button should be clear, not yellow.)*



**3. Slide the cap on until it stops.**

Align the notch on the cap with the notch on the device.



**4. Adjust the depth.**

Turn the Comfort Dial™ to start at 2 or 3. For tougher skin, dial to a higher number.



**5. Prime the lancet device.**

Press the plunger in as far as it will go, like a pen. The release button turns yellow when the device is ready.



**6. Hold the lancet device firmly against the side of your fingertip.**

Press the yellow release button. Remember, the lancet does not come out of the center of the cap.



**7. Advance to a new lancet.**

For virtually pain-free testing, always use a new lancet.

*(See other side for instructions on how to load a new lancet.)*

# Using your ACCU-CHEK® Multiclix Lancet Device

## Testing From a Site Other Than Your Fingertip

Some blood glucose meters test from sites other than your fingertip, such as palm, forearm, upper arm, thigh, or calf.\*

### 1. Slide the clear cap on by aligning the notch

on the cap with the notch on the lancet device. Be sure a lancet drum is loaded.

For palm testing, try a setting of 3, all other non-fingertip areas use a setting of 5½.

### 2. Prime the lancet device.

Press the plunger in as far as it will go, like a pen. The release button turns yellow when the device is ready.

### 3. Press the lancet device firmly against the skin.

Pump the lancet device slowly up and down to increase circulation.

### 4. Press yellow release button while maintaining firm, steady pressure to assist blood flow


until an adequate drop of blood forms.



1



2




Release button turns yellow

3



4



Press release button

## More Tips for Testing From a Site Other Than Your Fingertip

- Between steps 2 and 3, rub the skin prior to lancing to increase blood flow.
  - If you experience difficulty obtaining a sufficient drop of blood, contact our ACCU-CHEK Customer Care Service Center. We will guide you through the testing procedure step by step.
- \* Refer to your owner's booklet or contact your doctor to learn more before testing from sites other than your fingertip.

## For virtually pain-free testing, always use a new lancet.

### To load a new lancet:

**A.** Twist the plunger ¼ turn forward until it stops, then twist back all the way.

**B.** The lancet counter will decrease by one white bar. When no white bars are visible, you are using the final lancet.

### To remove the drum:

Take the cap off and pull the drum straight out. Once removed, a lancet drum cannot be reinserted.

### To insert a new drum:

See the other side for instructions.



Lancet Counter



The ACCU-CHEK Customer Care Service Center is available 24 hours a day at **1-800-858-8072**. For a virtual product demonstration, go to **accu-chek.com**