# Welcome to the TrueTrack Smart System® Blood Glucose Monitoring System

Congratulations on your purchase of the TrueTrack Smart System® Blood Glucose Monitoring System. You now have a very simple and accurate way to test your blood glucose (sugar) level, anytime, anywhere.

## Our Commitment to You

Our goal is to provide you with quality healthcare products and dedicated customer service. If you have questions about using TrueTrack Smart System® products, visit our web site at:  
**[www.homediagnostics.com](http://www.homediagnostics.com)**.

***Caution! Please read complete Owner's Booklet and all product Instructions for Use.***

<b>SYMBOLS:</b> Sterile	Control Level	<b>SN</b> Serial Number	Single Use Only	Number of Tests	Use By Date
For Assistance Call	Attention! Read Instructions for Use.	Storage Temperature Range	Authorized Representative		
Lot Number	For <i>in vitro</i> Diagnostic Testing Only	Code Number	Manufactured By		

## Importance of Blood Glucose Monitoring

The more you know about diabetes, the better you are able to take care of yourself. As you meet with your Doctor or Diabetes Healthcare Professional, you will work together to determine the target range for your blood glucose results and how often you should test. Having the majority of your results within your target range helps delay or prevent complications of diabetes.

Using the TrueTrack Smart System® Blood Glucose Monitoring System for self-testing shows you and your Doctor or Diabetes Healthcare Professional how well your daily treatment plan is working to control your blood glucose. Do not make changes to your treatment plan without discussing with your Doctor or Diabetes Healthcare Professional.

The TrueTrack Smart System® is intended for *in vitro* (outside the body) testing only, and is not for the diagnosis of diabetes or for testing blood glucose in newborns.

Use of the TrueTrack Smart System® Blood Glucose Monitoring System in a manner not specified in this Owner's Booklet is not recommended and may affect your ability to determine your blood glucose.

## Important Information

**Important:** For the best possible performance from TrueTrack Smart System<sup>®</sup> we recommend you:

- **Read all instructions** before you begin to test.
- **Perform Check Strip and Glucose Control Tests *before*** testing your blood.

*Note: Glucose Control may not always be included. Check Contents List on outside packaging. If Glucose Control is not listed, contact place of purchase.*

- **Set date and time** for correct Morning Average values.
- **Check Codes.** Code in Meter Display must match Code on vial of TrueTrack Smart System<sup>®</sup> Test Strips being used and Code Number printed on Code Chip.
- **Use only TrueTrack Smart System<sup>®</sup> Test Strips and TrueTrack Smart System<sup>®</sup> Glucose Control with TrueTrack Smart System<sup>®</sup> Meter.**
- **Use only fresh capillary whole blood samples** for testing. TrueTrack Smart System<sup>®</sup> is not recommended for use with venous samples. Meter displays results in ***Plasma*** values.
- **The TrueTrack Smart System<sup>®</sup> is recommended for testing of human whole blood only.**

## Health Related Information

- Severe dehydration and/or excessive water loss may cause false low results.<sup>2</sup> If you believe you are suffering from severe dehydration, contact your Doctor or Diabetes Healthcare Professional immediately.
- Test results equal or below 80 mg/dL (4.4 mmol/L) indicate low blood glucose (**hypoglycemia**).<sup>3</sup> Test results greater than 250 mg/dL (13.9 mmol/L), indicate high blood glucose (**hyperglycemia**).<sup>4</sup> Perform another blood test to assure that your results are correct. If you continue to get similar results and you are having symptoms that suggest glucose is too low or too high, follow the treatment advice given by your Doctor or Diabetes Healthcare Professional.

## **Table of Contents**

Fast Test Guide .....	inside front cover
Your Target Blood Glucose Range .....	inside back cover
Introductory Information .....	i-v
<b>Know Your System</b>	
Meter .....	1-3
Test Strip .....	4
Test Strip Vial Label .....	5
<b>Getting Started</b>	
Turning Meter On and Off.....	6
Meter Set Up (Time, Date, Units of Measure) .....	7-8
Coding .....	9-10
Quality Control Tests .....	11
Check Strip Test .....	11
Glucose Control Test.....	12-15
<b>Testing your Blood</b>	
Obtain a Blood Sample.....	16-17
Testing Blood Glucose.....	18-19
Result Warning Messages .....	20
<b>Meter Memory</b>	
Viewing Results in Memory .....	21
Morning Average (14 and 30 Day).....	22
<b>Care and Troubleshooting</b>	
Caring for TrueTrack Smart System® .....	23-24
Battery Installation .....	25
Troubleshooting.....	26
Display Messages .....	27-29
Meter Specifications .....	30
Chemical Composition.....	31
Meter Warranty .....	32
References.....	33

# Know Your System

## Meter


Top of Meter

Front of Meter


- ① **← Button** - Decrease numbers in Set Up; remove Control Symbol; move backward by date when viewing results in Memory.
- ② **S Button** - Select settings; view Morning Average values; view results in Memory.
- ③ **→ Button** - Increase numbers in Set Up; add Control Symbol; move forward by date when viewing results in Memory.
- ④ **Display** - Shows test results, messages and other information.
- ⑤ **Test Port** - Insert TrueTrack Smart System<sup>®</sup> Test Strip or Check Strip here.

**Right Side of Meter****Left Side of Meter****Back of Meter**

- ⑥ **Data Port** - Contact Customer Service at the number located on back cover.
  - ⑦ **Code Chip Port** - Insert the Code Chip here.
  - ⑧ **Battery Compartment** - Use one non-rechargeable 3V  Lithium Battery (#CR2032), positive (“+”) side up.
- Note: See Pg. 25 for details on discarding old battery.*
- ⑨ **Meter Label** - Contains telephone number for assistance.
  - ⑩ **Serial Number Label** - Identifies Meter when calling for assistance.

## Display

1. Time, Date, 14-30 Day Average, Code.
2. Test Result in Memory.
3. Morning Average value.
4. Check Strip Test correct.
5. Test results.


**Note:** *A decimal point appears when results are in mmol/L.*

6. Battery Symbol (see *Display Messages*).
7. Temperature Symbol (see *Display Messages*).
8. Apply Sample Symbol.
9. Control Symbol (see *Glucose Control Test*).
10. Glucose Units of Measure.
11. Time is P.M.

**Note:** *AM does not appear.*


## Test Strip


- ③ **Contact End** - Insert into Meter with Contacts (metallic blocks) facing up.
- ④ **Sample Tip** - Place sample (blood or Glucose Control) to edge of Tip.

## Placement of Sample to Test Strip

*Correct*

*Incorrect*

**Note:** *Do not apply blood or Glucose Control to top of Test Strip!*

*Use Strip quickly after removal from vial. Strips that have been left out of vial too long will give error message when used.*

## Test Strip Vial Label

<b>Code</b>	1027	①
<b>LOT</b>	EA645A4	
	28DEC07	②
0	89-121	
1	177-239	③
2	356-493	
mg/dL / mmol/L		
<b>May 13, 2007</b>		②


- ① **Code** - Match Code in Meter Display with Code on Code Chip *and* Code on Test Strip vial label *before* testing. All Codes **MUST** match for accurate results.
- ② **Use by Dates** - Write date first opened on vial label. Discard vial and unused Strips 120 days after written date or date printed next to  symbol, whichever comes first.
- ③ **Control Level Range** - Range of numbers in which Glucose Control Test result must fall to assure System is working properly.

**Caution!** *Use of Test Strips or Glucose Control past the Use by Dates may give inaccurate or incorrect test results. Discard out-of-date products and test with new products.*

# Getting Started


## To Turn Meter On...

- ① For Set Up, press and hold **S** until full display appears.
- ② For Check Strip Test, insert Check Strip into Test Port. For Blood or Control Test, insert Test Strip into Test Port.
- ③ For Morning Average/Memory, press and release **S**.

## To Turn Meter Off...

- ① After viewing Set Up or Morning Average/Memory, press and hold **S** until Display is blank.
- ② After Testing, remove Test Strip or Check Strip from Test Port.
- ③ Meter turns off after two minutes of non-use.

# Meter Set Up


1. Start with Meter off. Do not insert Test Strip.

Press and hold **S**. Date/Time appears in display. Release **after** full display appears (*shown at right*) and Meter beeps.

## Set Hour / Minutes


2. The hour flashes.  
Press **➤** to increase or **➤** to decrease. Press **S** when correct.

**Note:** To set “PM”, scroll until “PM” is below correct time. “PM” is not an option if factory setting is for a 24-hour clock.


3. The minutes flashes.  
Press **➤** to increase or **➤** to decrease. Press **S** when correct.

## Set Month/Day/Year


- The month (number) flashes. Press **+** to increase or **-** to decrease. Press **S** when correct.
- The day (number) flashes. Press **+** to increase or **-** to decrease. Press **S** when correct.
- The year flashes. Press **+** to increase or **-** to decrease. Press **S** when correct.

## Set Unit of Measure

- The unit of measure (mg/dL or mmol/L) flashes. Press **+** or **-** to change. Press **S** to select.

**Note:** *mg/dL is default setting for Meters sold in United States.*

OR

## Exit Set-up Mode

To review settings, press **S**. Make changes as needed. Press and hold **S** until Display goes blank. Options are saved. Meter turns off.

## Coding

Code Chip is packaged with each box of Test Strips.

**Note:** *After all Test Strips are used from box, remove Code Chip from Meter and discard. Open new box of Test Strips and insert new Code Chip into Meter before testing.*

1. Start with Meter off. Remove and open Code Chip package from box of Test Strips to be used. Match Code printed on Code Chip to Code printed on vial label.

**Caution!** *If Codes do not match, do not use Test Strips. Call for assistance.*


2. If Codes match, insert Code Chip into Code Chip Port as shown, locking in place.

3. Insert Test Strip into Test Port.  
Meter turns on. “Code” and number appear in Display.

Code 1027

4. Match to Code number on vial label. If Codes match, begin testing.


**Caution!** If “Code ----” shows in Display, reinsert Code Chip and then reinsert Test Strip. If “Code ----” still shows in Display or if Codes do not match, **DO NOT** use Test Strips. Test results could be incorrect. Call for assistance.

# Quality Control Tests

## Check Strip Test

Check Strip is stored in Carrying Case. Check Strip Test is used to check that Meter is working properly. Perform a Check Strip Test when:

- Using Meter for first time.
  - Results do not match how you feel.
  - Meter was dropped or banged.
1. Start with Meter off. Remove Check Strip from Carrying Case. Insert Check Strip Contact End ① into Test Port with dimple facing up. Meter turns on.
  2. Check Display to make sure that all segments are showing. (See Pg. 3.)
  3. After a few seconds, Meter displays OK ✓ if test is correct.

**Caution!** *If segments are missing from Display or OK ✓ does not appear, do not use Meter for testing. Call for assistance.*

4. Remove Check Strip. Return to Carrying Case. Meter turns off.

## Glucose Control Test

Glucose Control Test uses TrueTrack Smart System<sup>®</sup> Glucose Control to check performance of System.

**Note:** *Glucose Control may not always be included. Check Contents List on outside packaging. If Glucose Control is not listed, contact place of purchase or Technical Service at the number located on back cover to obtain Glucose Control.*

Perform Control Tests when...

- Using Meter for first time. We recommend you perform several Control Tests for practice before testing with blood.
- Opening a new vial of Test Strips.
- Meter results seem unusually low or high.
- Meter was dropped.

**Caution!** *Result ranges are printed on Test Strip vial label. These ranges are for Control Test results only and **are not** recommended result ranges for your blood glucose.*

***Do not ingest Glucose Control.***

Most current Control Test result is stored in Memory.


## How to Test Glucose Control

1. Wash hands with soap and water, dry thoroughly. Allow Control, vial of Test Strips and Meter to adjust to room temperature (68-77°F/20-25°C). Write date first opened on Control Label.
2. Check dates on Control label and Test Strip vial label. Do not use Control if 3 months past written opened date *or* past date printed next to  symbol. Do not use Test Strips 120 days past written opened date or past date printed next to  symbol. Discard out-of-date products and use new products if either date has passed.


3. Swirl or invert bottle gently to mix Control. **DO NOT SHAKE!**
4. Remove Test Strip from vial. Recap vial immediately. Wipe top of vial cap with a wet tissue, dry thoroughly.

**Note:** Use Strip quickly after removal from vial. Strips that have been left out of vial too long will give an error message. Discard old Strip and use new Strip for testing.

5. Insert Test Strip into Test Port. Meter turns on.

6. Match Code in Display with Code on Test Strip vial label being used. If Codes do not match, see *Coding* (pg. 9).
7. Press . Control Symbol appears in Meter Display.

**Note:** If you decide not to perform a Control Test, press  to remove Control Symbol.

8. Turn Control bottle upside down. Squeeze one drop of Control onto a clean tissue. Wipe off bottle tip.

9. Gently squeeze a drop onto top of *closed* vial cap. Touch edge of Test Strip Sample Tip to drop of Control (See pg. 5).

**Note:** A small piece of unused aluminum foil or clear plastic wrap may be used instead of vial cap. Dispose after use.

10. Allow drop to be drawn into Strip. Remove Strip from drop when Meter beeps.


**Note:** *If Meter does not beep or begin countdown soon after placing Sample Tip to Control drop, discard Strip. Repeat test with new Test Strip. If problem persists, see Troubleshooting.*

- 11.** Compare result to Control Level range printed on Test Strip vial label:

If result is in range, System can be used for testing blood. If result does not fall within range, repeat test using a new Strip.

0	89-121
1	177-239
2	356-493
mg/dL / mmol/L	

**Caution!** *If Control Test result continues to be outside range, System should not be used for testing blood. Call for assistance.*

- 12.** Remove Test Strip from Meter and discard. Meter turns off.
- 13.** Wipe Control drop off of vial cap with a clean, wet tissue and dry thoroughly.

## How to Obtain a Blood Sample

Refer to Lancing Device “Instructions for Use” for detailed instructions.

### From Fingertip


1. Prepare fingertip by washing hands in warm, soapy water. Rinse well. Dry thoroughly.
2. Place end of Lancing Device against tip of finger. Lance fingertip.
3. Set Lancing Device aside. To help blood drop form, lower hand to waist level, gently massaging finger from palm to fingertip. Allow blood drop to form before attempting to apply to Test Strip.

Always remove and discard used Lancet in appropriate container when testing is complete.

## For Forearm Testing

1. Select area to be lanced. Wash with soap and warm water, rinse and dry thoroughly.
2. Rub area vigorously or apply a warm compress to stimulate blood flow.
3. Place tip firmly against forearm. Press trigger button. Apply firm pressure on lancing device until an adequate blood drop is formed.

### *Important Notes Regarding Forearm Samples Used in Testing Blood Glucose:*


- Check with your Doctor or Diabetes Healthcare Professional to see if forearm testing is right for you.
- Blood glucose results from forearm are not always the same as results from fingertip. Do not use forearm for testing when blood sugar is rapidly rising or falling, e.g., within 2 hours of eating, after taking insulin, immediately after exercise, or when you are ill or under stress. Use fingertip sample for testing.
- Do not use forearm samples for testing if you think your blood sugar is low, if you don't regularly experience symptoms when your blood glucose is low, or if your results from forearm samples do not match the way you feel. Use fingertip sample for testing.<sup>5</sup>

## How to Test Blood Glucose

1. Check dates on Test Strip vial being used. Do not use if 120 day past written date or date next to  symbol, whichever comes first
2. Wash hands. Rinse well and dry thoroughly.
3. Remove Test Strip from vial. Recap vial immediately.

**Note:** Use Test Strips quickly after removal from vial. Strips that have been left out of vial too long will give an error message. Discard old Strip and use new Strip for testing.

4. With Meter off, insert Test Strip Contact End (Contacts facing up) into Test Port. Meter turns on.
5. Match Code in Display to Code on Test Strip vial being used. If Codes do not match, see *Coding*.
6. Lance fingertip or forearm. Allow drop to form.


7. Touch edge of Test Strip Sample Tip to blood drop and hold. (See pg. 5.) Allow blood to be drawn into Test Strip until Meter beeps.

**Note:** *If Meter did not beep or begin countdown soon after touching blood drop to Test Strip Sample Tip, discard Strip. Repeat test with new Test Strip and new blood drop. If problem persists, see Troubleshooting.*

8. After countdown, result is displayed with date and time. Remove Test Strip and discard. Meter turns off.

Result is stored in Meter's Memory.  
Record result with date in log book.

## Result Warning Messages

### WARNING!

If blood glucose test result is less than 20 mg/dL (1.1 mmol/L), “**LO**” appears in Meter Display.

If blood glucose test result is greater than 600 mg/dL (33.3 mmol/L), “**HI**” appears in Meter Display.

Always repeat test to confirm Low (“**LO**”) and High (“**HI**”) results. If results still display “**LO**” or “**HI**”, call your Doctor or Diabetes Healthcare Professional *immediately*.

**Note:** “**LO**” results are included in the Morning Average as 20 mg/dL (1.1 mmol/L). “**HI**” results are included as 600 mg/dL (33.3 mmol/L).

## Viewing Results in Memory


Memory stores 365 results, which are displayed from most recent to oldest.

Start with Meter off. Press and release **S**. Meter displays date and time followed by Morning Average. Press **S** again to view most recent result in Memory.

To scroll results forward, press **➤**. To scroll results backward, press **➤**.

Latest Glucose Control Test result appears after earliest blood test result and shows Control Symbol in Display. Only one Glucose Control Test result is stored in Memory.


## Morning Averages (14 and 30 Day)

For Morning Average values, at least one blood test must be performed between 4:00 am - 9:59 am according to Meter clock. Meter clock must be set correctly (See *Set-Up*) for accurate Morning Average values.

Start with Meter off. Press and release **S**. Display shows 14 day Average value and then 30 day Average value.

Meter switches between 14 and 30 day Average values for 2 minutes before turning off.

If there are no Morning Average values, three dashes are displayed.

## Caring for TrueTrack Smart System®

- Store System (Meter, Check Strip, Glucose Control, Test Strips, Code Chip) in Carrying Case to protect from liquids, dust and dirt.
- Store System in a dry place at room temperature (36°-86°F / 2°-30°C). **DO NOT FREEZE.**

### Meter Care

- Wipe Meter with clean, lint-free cloth dampened with one of the following:
  - Mild detergent or mild soap and water
  - 10% household bleach and water
- Never immerse Meter in liquids or allow any liquids into Test and Data Ports.

### Check Strip Care

- Do not bend Check Strip.
- Do not apply blood, Glucose Control or any other liquid to Check Strip.
- To clean, wipe with damp cloth, dry thoroughly.
- Do not write on Check Strip.

### Glucose Control Care

- Write date opened on Control label. Discard 3 months after opening or after date printed next to  symbol, whichever comes first.
- After use, wipe bottle tip clean and recap tightly.
- Store at room temperature (36°-86°F / 2°-30°C).
- Swirl bottle to mix. Do not shake.

## Test Strip Care


- Store Strips in original vial only.
- Write date opened on Strip vial. Discard unused Strips from vial 120 days after opening or after date printed next to  symbol, whichever comes first. Use of Strips past either date may give incorrect results.
- Close vial immediately after removing Strip. Use Strip quickly after removal from vial. Never store Strips outside of vial. Store in a dry place at room temperature below 86°F/30°C. **DO NOT FREEZE.**
- Do not reuse Strip.
- Apply only fresh capillary blood or TrueTrack<sup>®</sup> Glucose Control to Strips.
- Do not bend, cut or alter Strips in any way.

## Code Chip Care

- Code Chip must remain in Meter until all Strips from box are used, then discard old Code Chip.
- Insert new Code Chip from new box of Strips immediately before using Test Strips.

## Battery Installation

Replace Battery when Low or Dead Battery Symbols appear in Display or Meter does not turn on.

1. Lift tab on Battery Cover.
2. Turn Meter over, tap gently to loosen and remove Battery.
3. Insert new Battery, positive (“+”) side up. Close cover.

**Note:** Use non-rechargeable 3V  Lithium Battery (#CR2032).

4. Discard old battery in appropriate container.
5. Turn Meter on. If Meter will not turn on, check that Battery was installed properly. If not, reinsert Battery and try again. Call for assistance if problem persists.

**Caution!** Batteries might explode if mishandled or incorrectly replaced. Do not dispose of battery in fire. Do not disassemble or attempt to recharge Battery. Dispose according to local/country specific regulations.

**Note:** Replacing battery may affect date and time settings. Check date and time by pressing  and releasing. To reset, see Set Up.

# Troubleshooting

## *1) After inserting Test Strip, Meter does not turn on.*


<b>Reason</b>	<b>Action</b>
Strip inserted upside down or backwards	Remove Strip. Re-insert correctly.
Strip not fully inserted	Remove Strip. Re-insert Strip fully into Meter.
Strip Error	Repeat with new Strip.
Dead or no battery	Replace battery.
Battery in backwards	Battery positive (“+”) side must face up.
Meter Error	Call for assistance.

## *2) After applying sample, test does not start/ Meter does not beep or begin countdown.*

<b>Reason</b>	<b>Action</b>
Sample drop too small	Repeat test with new Strip and larger drop.
Sample applied after two minute shut-off	Repeat test with new Strip and apply drop after Meter beeps.
Problem with Strip	Repeat with new Strip.
Problem with Meter	Call for assistance.

## Messages

### Display


### Reason

**Temperature Error**  
Temperature change too quick


### Action

Wait 10 minutes for Meter and Strips to reach room temperature before testing.

**Temperature Too Cold**  
Meter temperature less than 50°F (10°C)

Move Meter and Test Strips to area between 50°-104° (10°- 40°C) before testing.


**Temperature Too Hot**  
Meter temperature greater than 104°F (40°C)

Move Meter and Test Strips to area between 50°-104° (10°- 40°C) before testing.

**Sample Not Detected**

Retest with new Test Strip and larger sample.



**Used Test Strip or Test Strip outside of vial too long.**




Repeat with new Test Strip. If error persists, call for assistance.

## Messages *(continued)*

<u>Display</u>	<u>Reason</u>	<u>Action</u>
	<b>Meter Error</b>	Call for assistance.
	<b>Test Strip Error</b>	Retest with new Test Strip. If error persists, call for assistance.
	<b>Strip Removed During Test</b>	Retest with new Strip. If error persists, call for assistance.
	<b>Check Strip Error</b>	Repeat test. If error persists DO NOT use Meter for testing. Call for assistance.
	<b>Code Chip Error</b>	Remove Code Chip and reinsert. If error persists, call for assistance.

## Messages *(continued)*

<u>Display</u>	<u>Reason</u>	<u>Action</u>
	<b>Communication Error</b>	Turn Meter off and on again. If error persists, call for assistance.
	<b>Low Battery</b>	About 50 tests can be done before Battery must be replaced.
	<b>Dead Battery</b>	Replace Battery.
	<b>Code Chip Error</b>	Remove Code Chip and re-insert. If error persists, call for assistance.
 	<ul style="list-style-type: none"><li>• <b>WARNING!!</b> <b>HI - High Results</b> Result more than 600 mg/dL (33.3 mmol/L)</li><li>• <b>LO - Low Results</b> Result less than 20 mg/dL (1.11 mmol/L)</li></ul>	<ul style="list-style-type: none"><li>• <b>WARNING!!</b> Retest with new Strip. If result is still “HI” or “LO” contact Doctor <i>immediately.</i></li></ul>

## Meter Specifications

**Result Range:** 1.1 - 33.3 mmol/L (20 - 600 mg/dL)

**Sample Size:** Minimum 1 microliter (1  $\mu$ l)

**Sample:** Fresh capillary whole blood or Glucose Control

**Test Time:** 10 seconds

**Result Value:** Plasma values

**Assay Method:** Electrochemical

**Power Supply:** One 3V  Lithium Battery  
#CR2032 (non-rechargeable)

Total power when active  
at full battery = 8.6mW

**Battery Life:** Approximately 1,100 tests or 1 year

**Automatic shut-off:** After two minutes of non-use

**Weight:** 47 grams (1.66 ounces)

**Size:** 89 x 55 x 17mm (3.52" x 2.15" x 0.67")

**Memory Size:** 365 blood glucose results,  
1 Control Test result

## **System Operating Range (Meter & Test Strips):**

**Relative Humidity:** 10-90% (Non-condensing)

**Temperature:** 10°- 40°C (50°-104°F)

**Hematocrit:** 30-55%

**Note:** *Use within specified environmental conditions only.*

## Chemical Composition

**Test Strips:** Glucose Oxidase (*Aspergillus niger*)  
2.5 Units, Mediators, Buffers and Stabilizers.

**Glucose Control:** Contents: Volume: 3 mL. Water - 73%, D-glucose - 0.09 - 0.20%, Viscosity Enhancing Agent - 25%, Inorganic Salts - 1.8%, Amaranth - 0.08%, and Preservatives - 0.03%.

## **Limited Warranty**

Home Diagnostics, Inc. provides the following Warranty to the original retail purchaser of the TrueTrack Smart System® Meter:

- 1) Home Diagnostics Inc. warrants this Meter to be free of defects in materials and workmanship at the time of purchase. If the meter is inoperative prior to five (5) years from the date of the original retail purchase, Home Diagnostics, Inc will replace the Meter with an equivalent Meter, at its option, at no cost to the purchaser. Failure of the meter due to abuse or use not in accordance with the instructions for use is not covered by this Warranty.
- 2) This Warranty does not include the battery supplied with the Meter.
- 3) Do not take the Meter apart. This action will void the Warranty and cause the Meter to display false results.
- 4) The duration of any implied Warranty, including any implied Warranty of merchantability or fitness for a particular purpose shall be limited to five (5) years from the date of original retail purchase. Some states do not allow limitations on how long an implied Warranty last, so the above limitation may not apply.
- 5) Home Diagnostics, Inc. disclaims liability for incidental or consequential damages for breach of any expressed or implied Warranty, including any implied Warranty of merchantability or fitness for a particular use with respect to the Meter. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations or exclusion may not apply.
- 6) This Warranty gives the user specific legal rights, and the user may also have other rights which vary state to state.

Your TrueTrack Smart System® Technical Service Representative will be able to provide detailed information regarding procedures for returning your Meter, if necessary.

## References

1. Joslin Diabetes Center & Joslin Clinic, Inc., (Sept. 2001). *Guidelines for Pharmacological Management of Type 2 Diabetes*. [Electronic Version]. Retrieved August 22, 2002 from [www.joslin.harvard.edu/education/library/wbgoal.shtml](http://www.joslin.harvard.edu/education/library/wbgoal.shtml).
2. Hackett, J., Aziz, K., *Review Criteria Assessment of Portable Blood Glucose Monitoring in vitro Diagnostic Devices Using Glucose Oxidase, Dehydrogenase or Hexokinase Methodology*; FDA Draft Document, p. 6 (1997).
3. Joslin Diabetes Center, (2002). *What Can I Do to Prevent Serious Hypoglycemic Episodes When I Am Hypoglycemic Unaware?* [Electronic Version]. Retrieved August 22, 2002 from [www.joslin.org/education/library/hypoglycemic\\_unaware.shtml](http://www.joslin.org/education/library/hypoglycemic_unaware.shtml).
4. Joslin Diabetes Center & Beth Israel Deaconess Medical Center, (2002). *Guidelines for Management of Hypoglycemic Emergencies* [Electronic Version]. Retrieved August 22, 2002 from [www.diabetesmanagement.joslin.org/-guidelines/hyperglycemic-guide.pdf](http://www.diabetesmanagement.joslin.org/-guidelines/hyperglycemic-guide.pdf).
5. Food and Drug Administration 2003. *Diabetes Information - Glucose Meters and Diabetes Management*. [Electronic Version]. Retrieved August 11, 2003 from [www.fda.gov/diabetes/glucose.html](http://www.fda.gov/diabetes/glucose.html).

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