AW 06488101B

**ONETOUCH<sup>®</sup>**  
**Select<sup>™</sup>**  
**Control Solution**

### Intended Use

For use with the OneTouch<sup>®</sup> Select<sup>™</sup> Blood Glucose Meter to **check that the meter and test strips are working together properly** and that you are performing the test correctly.

### Summary and Explanation

OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution contains a known amount of glucose that reacts with OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips. The test result you obtain should fall within the range of results that you will find printed on the vial of the test strips that you are using.

**⚠ CAUTION:** The control solution range printed on the test strip vial is for control solution only. It is used to test meter and test strip performance. **It is not a recommended range for your blood glucose level.**

**⚠ WARNING:** Keep the meter and testing supplies away from young children. Small items such as the battery door, batteries, test strips, lancets, protective disks on the lancets, and control solution vial cap are choking hazards.

### Do a control solution test:

- To practice the test process, instead of using blood.
- Once a week.
- Whenever you open a new vial of test strips.
- If you suspect that the meter or test strips are not working properly.
- If you have had repeated unexpected blood glucose results.
- If you drop or damage the meter.

A single control solution test that is within the expected range is enough to tell you that the meter and test strips are working together and that you are doing the test properly.

### Composition

The control solution consists of glucose in water with buffers and stabilizers, a viscosity adjusting agent, a preservative, and a dye. The approximate glucose concentration for the OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution is 0.11%. An expected range of glucose results for the control solution is printed on the label of your OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips.

### Important Safety Information

- Use only OneTouch<sup>®</sup> Select<sup>™</sup> Control Solutions with the OneTouch<sup>®</sup> Select<sup>™</sup> Meter. Write the discard date (date opened plus three months) in the space provided on the vial label.
- **Do Not** use after the expiration or discard date, whichever comes first.
- For *in vitro* diagnostic use.

**⚠ CAUTION:** Do Not swallow control solution; it is not for human consumption. Do Not apply control solution to the skin or eyes as it may cause irritation.

### Storage and Handling

Close the bottle tightly and store at room temperature below 86°F. **Do Not** refrigerate.

### Test Procedure

Read the instructions carefully in your owner's booklet to perform the control solution test.

Materials provided: OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution

Materials required but not provided: your meter, owner's booklet, OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips

Follow these steps whenever you test with OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution:

1. Control solution tests must be done at room temperature (68–77°F). Make sure your meter, test strips, and control solution are at room temperature before testing.
2. Shake the control solution vial well before testing.
3. Discard the first drop of control solution you squeeze out and wipe the dispenser tip.
4. Be sure that APPLY CONTROL appears on the meter display before applying control solution to the test strip. (See your owner's booklet for details.)
5. Apply the drop of control solution to the narrow channel in the top edge of the test strip.
6. Your OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips have two control ranges printed on the vial. Be sure to compare the result to the correct control range.

Tightly close the cap on the control solution vial immediately after use to avoid contamination.

Each vial contains sufficient volume for a minimum of 75 tests.

### Expected Results

At room temperature, OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution test results should fall within the expected range printed on the test strip vial label about 95% of the time.

If test results fall outside the expected range, first repeat the control solution test. Results outside the expected range may indicate:

- An error in the way you are doing the test.
- Expired or contaminated control solution.
- Damaged or expired test strips.
- Use of a test strip or control solution past its discard date.
- A problem with the meter.
- Control solution that is too warm or too cold.
- Failure to shake the control solution vial well.
- Failure to discard the first drop of control solution and to wipe the dispenser tip.
- Failure to match the result to the correct control range.

Refer to your owner's booklet for assistance. If you are unable to resolve the problem, contact Customer Service.

Used test strips may be considered biohazardous waste in your area. Be sure to follow your local regulations for proper disposal.

**⚠ CAUTION:** If you continue to get control solution test results that fall outside the range printed on the test strip vial, **Do Not** use the meter, test strips, or control solution. Call LifeScan Customer Service at 1 800 227-8862.

### OUR COMMITMENT TO YOU

Our goal is to provide you with quality health care products and dedicated customer service. If you are not fully satisfied with this product, or if you have questions about the use of any LifeScan product, please call us at 1 800 227-8862.

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